

Whole-Mount Embryo labeling with (Monoclonal) Antibodies

Source: Zoltan Varga; this protocol is a derivative of the chapter found in the Zebrafish Book, Westerfield, 2000/ Karen Larison)

Solutions

PO₄ buffer (0.1 M, pH 7.3): 80 ml 0.1 M Na ₂ HPO ₄ (13.8 g Na ₂ HPO ₄ · H ₂ O/Liter ddH ₂ O) 20 ml 0.1 M NaH ₂ PO ₄ (26.8 g NaH ₂ PO ₄ · 7H ₂ O/Liter ddH ₂ O)	2X PBS: 8.0 g NaCl 0.2 g KCl 200 ml 0.1 M PO ₄ Buffer, pH 7.3 300 ml dH ₂ O When diluted 1:1 with dH ₂ O, final concentrations are 0.8% NaCl, 0.02% KCl and 0.02 M PO ₄ .	4% Paraformaldehyde (Fixative): Mix equal amounts of 8% paraformaldehyde and 2x Phosphate buffered saline (PBS, pH 7,3)
PBS/BSA/Triton: 50 ml 2X PBS, pH 7.3 1 g BSA 1 ml 10% Triton X-100 Add dH ₂ O to a total volume of 100 ml.	Block: 2% Normal goat serum (NGS) in PBS/BSA/Triton Primary antibody: 2% Normal goat serum (NGS) in PBS/BSA/Triton + 1° Antibody (dilute according to recommendation) Secondary antibody: (fluorescence, or enzyme tagged) 2% Normal goat serum (NGS) in PBS/BSA/Triton + 2° Antibody (dilute according to manufacturer's recommendations)	

Embryos remain in Eppendorf tubes throughout the staining procedure; remove and add solutions as indicated.

1. Fix embryos in 4% paraformaldehyde for 3-4 hours at room temperature or overnight at 4°C.
2. Wash 3x in 0.1 M PO₄ buffer, pH 7.4, for 5 min.
Permeabilization steps 3. – 6. are optional:
 3. Wash in dH₂O for 5 min.
 4. Freeze in acetone at -20°C for 7 min. to permeabilize tissue.
 5. Wash in dH₂O for 5 min.
 6. Wash in 0.1 M PO₄ buffer for 5 min.
7. Treat with PBS/BSA/Triton + 2% NGS for 30 min to block non-specific binding sites.
8. Soak in primary antibody (Mouse IgG) diluted in PBS/BSA/Triton + 2% NGS overnight at 4°C. Best in cold room on shaker to facilitate penetration. Shorter incubation times are also possible at room temperature.
9. Wash 4x 30 min in PBS/BSA/Triton at room-temperature.
10. Soak in secondary antibody diluted in PBS/BSA/Triton + 2% NGS overnight at 4°C or 5 hr at RT.
11. Wash 4x 30 min PBS/BSA/Triton (overnight at 4°C).
12. Transfer embryos into PBS.
13. Mount tissue in 4% Methyl cellulose on slide, using pieces of coverslips to serve as spacers.
The embryos can also be viewed using depression slides.

Thorough washes/rinses seem to be the most important factor to reduce background (if any) additional or extended wash/rinsing steps can always be included.

Instead of Acetone you can use **ProK** to permeabilize tissue (optional):

14. use PBS buffer, 10 $\mu\text{g/ml}$ (working dilution) ProK in ddH₂O. The exposure time will vary depending on the embryonic larval stage, the batch of Prok you use, and whether you have frozen and thawed the ProK stock solution several times. It is recommended to test each batch of ProK on embryos/larvae that you won't need for experiments, then aliquot and freeze your ProK stock at $-20\text{ }^{\circ}\text{C}$.

After ProK treatment

15. rinse at least 3x rapidly in ddH₂O,
16. rinse at least 3x 5min in PBS,
17. re-fix samples in 4% PFA for 20 min at RT,
18. rinse 3x 5min in PBS.