ZEBRAFISH INTERNATIONAL RESOURCE CENTER

Last Updated March 2, 2008

kita^{b5}

Nature of the mutation

The b5 allele contains a single nucleotide deletion (a deletion of T) at codon 846 of the kita gene, resulting in a frameshift and a premature stop codon (Parichy et al., Development 126: 3425-3436, 1999).

Genotyping assay

Genotyping of the b5 allele is based on the RFLP assay (Restriction Fragment Length Polymorphism; Botstein et al., Am. J. Hum. Genet. 32: 314-331, 1980). This method is used to detect a mutation that either creates or abolishes a site recognized by a specific restriction enzyme. In the RFLP assay, a sequence of interest is first PCR-amplified and then the PCR product is subjected to restriction enzyme digestion. The presence or absence of the mutation is determined by the resulting restriction pattern. The b5 mutation abolishes a site recognized by the Tsp509I restriction enzyme.

Primers:

KIT 03: 5' AGA CAA CCA CAA GCA CTT TAC CC 3' KIT 04: 5' CCA TAG ATG TTG ACC CAG ACT GG 3'

PCR program (55_30_30):

- 1. 94°C for 3 min
- 2. 94°C for 30 sec
- 55°C for 30 sec 3.
- 4 72°C for **30** sec
- 5. Go to step 2 (above) for 39 cycles
- 72°C for 5 min 6.
- 7. 8.0°C hold
- 8. **END**

Product size: 269 bp

Digestion of the PCR product with the Tsp509I restriction enzyme:

Product type	Product digestion	DNA fragments after digestion (bp)
PCR product derived from the WT template	cleaved	202 bp and 67 bp
PCR product containing the mutation	unaffected	269 bp

IMPORTANT NOTE: It is highly recommended to use WT positive controls to monitor whether enzyme digestion has been carried out to completion. Without this control, partially digested WT samples can be mistakenly regarded as heterozygous samples.